

**ELECTRONIC SOMATIC CELL COUNT**

**Fossomatic™ Minor  
(Raw Cow Milk, Raw Sheep Milk, Raw Goat Milk and Raw Water Buffalo Milk)  
IMS #16**

**(Unless otherwise stated all tolerances ±5%)**

- 1. Laboratory Requirements (see Cultural Procedures (CP) items 33 & 34)** \_\_\_\_\_
- a. Un-preserved samples may be tested up to 72 hours after initial collection \_\_\_\_\_
- b. Samples may be tested up to 7 days after initial collection if preserved with 0.02% 2-bromo-2-nitropropane-1,3-diol (Bronopol™) or 0.05% potassium dichromate (K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub>) \_\_\_\_\_

**PRE-REQUISITE**

- 2. Comparative Test with DMSCC (co-requisite for certification)** \_\_\_\_\_
- a. Analyst(s) certified for DMSCC \_\_\_\_\_
- b. Each analyst seeking certification for the ESCC test shall perform the comparative test \_\_\_\_\_
  - 1. Test 4 samples (100K-200K, 300K-500K, 600K-800K and 900K-1.2M) in triplicate for both DMSCC (three separate smears each) and ESCC \_\_\_\_\_
  - 2. Results must be evaluated by State/Federal LEO and shown to be acceptable prior to official use of test in laboratory \_\_\_\_\_
  - 3. Copy of comparison and results in QC record (or easily accessible on file in the laboratory); kept for as long as analyst is certified \_\_\_\_\_

**APPARATUS**

- 3. See CP items 1-4** \_\_\_\_\_
- 4. Water Bath** \_\_\_\_\_
  - a. Circulating and thermostatically controlled to 37-42°C \_\_\_\_\_

**REAGENTS**

- 5. Reagents**
  - a. Dye Solution      Lot #: \_\_\_\_\_      Exp. Date: \_\_\_\_\_ \_\_\_\_\_
  - b. Clean 1            Lot #: \_\_\_\_\_      Exp. Date: \_\_\_\_\_ \_\_\_\_\_
  - c. Clean 2            Lot #: \_\_\_\_\_      Exp. Date: \_\_\_\_\_ \_\_\_\_\_

**6. Preparation**

- a. Ready to Use Dye Solution: Pour into a clean glass container designated for the Dye Solution (item 5.a). Use within 4 weeks of dispensing into container

Date Dispensed: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

- b. Clean 1 Solution: In a clean glass container, mix one unit (20 mL) Clean 1 (item 5.b) with deionized (DI) or Microbiologically Suitable (MS) water to make 1 liter, store and use within 4 weeks; when stored at 2-8°C, use within 8 weeks

Date Prep.: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

- c. Ready to Use Clean 2 Solution: Pour into a clean glass container designated for the Clean 2 Solution (item 5.c). Use within 4 weeks of dispensing into container

Date Dispensed: \_\_\_\_\_ Exp. Date: \_\_\_\_\_

**7. All solutions labeled with date prepared and expiration date**

**START UP**

**8. Cell Count**

- a. Check that the volume of Dye, Clean 1 and Clean 2 solutions in the supply containers is sufficient for the number of samples to be tested
- b. Solutions not used beyond expiration date(s)
- c. Perform the "Start Up" Job sequence: If the Zero Count is > 6, repeat "Clean Cuvette" and re-check the zero

- d. **IF ANY ABOVE PARAMETERS ARE OUT OF VARIANCE, CORRECT BEFORE PROCEEDING**

- e. Maintain records on all parameters each time instrument is used

**9. Milk Standards**

- a. Commercially prepared: \_\_\_\_\_

Lot #: \_\_\_\_\_ Date Rcd: \_\_\_\_\_

1. Four standards in ranges 100K-200K, 300K-500K, 600K-800K and 900K-1.2M
2. Perform DMSCC in triplicate on each standard in set and average counts; maintain records
3. Perform DMSCC check in rotation by all certified analysts

4. Use standards within one week \_\_\_\_\_

Lab Exp. Date: \_\_\_\_\_

b. Certified provider: \_\_\_\_\_

Lot #: \_\_\_\_\_ Exp. Date: \_\_\_\_\_ Date Rcd: \_\_\_\_\_

1. Four standards in ranges 100K-200K, 300K-500K, 600K-800K and 900K-1.2M \_\_\_\_\_

2. Maintain copies of all provided DMSCC values \_\_\_\_\_

3. Measure and maintain records of temperature (0.0-7.5°C) of standards as received \_\_\_\_\_

4. Maintain copies of all correspondence regarding problems \_\_\_\_\_

5. Standards used by manufacturer's expiration date \_\_\_\_\_

c. Laboratory prepared (weekly) \_\_\_\_\_

1. Prepare from raw milk > 18 hours old preserved with 0.05% potassium dichromate ( $K_2Cr_2O_7$ ) \_\_\_\_\_

2. Or, preserved with 0.02% 2-bromo-2-nitropropane-1,3-diol (Bronopol™) \_\_\_\_\_

3. Standards cannot be preserved with formalin \_\_\_\_\_

4. Prepare 4 standards in ranges 100K-200K, 300K-500K, 600K-800K and 900K-1.2M; use within one week \_\_\_\_\_

Lab Prep. Date: \_\_\_\_\_ Lab Exp. Date: \_\_\_\_\_

5. Perform DMSCC in triplicate on each standard and average counts; maintain records \_\_\_\_\_

6. Perform DMSCC check in rotation by all certified analysts \_\_\_\_\_

d. Hourly Control Sample (instrument drift check) \_\_\_\_\_

1. Use one of the standards (items 9.a, b or c) in the 600-800K range, test in triplicate and determine average \_\_\_\_\_

2. Optionally, prepare sufficient control/sample 600-800K range, test in triplicate and determine average \_\_\_\_\_

## PROCEDURE

### 10. Testing Standards (each time instrument used)

- a. Heat standards to 37-42°C (using a temperature control) and test within 30 min of reaching temperature, use once and then discard, i.e. do not re-use
- b. Mix by inverting at least 10x, test standards within 3 min
- c. Test the standards in triplicate and average the counts for each level; maintain records
- d. Each standard's average must be within 10% of the DMSCC (item 9) for that level, except within 15% for 100K-200K standard; maintain records
- e. Repeatability - a standard in the 300K to 800K range must have a coefficient of variation ( $C_V$ ) of 5% or less on 10 replicates (**Refer to Operating Manual**); maintain records
- f. Alternatively, set and run standard check as a "Custom Job", enter DMSCC values (item 9) into Excel™ macro, starting the job will enable 10.c through 10.e to be run and calculated automatically
- g. **THESE PARAMETERS MUST BE ACHIEVED BEFORE PROCEEDING**

### 11. Testing Samples

- a. Heat samples to 37-42°C (using a temperature control) and test within 30 min of reaching temperature
- b. Test samples within 10 min after removal from water bath
- c. Mix by inverting at least 10x, test samples within 3 min
- d. Record number of cells counted for each sample

### 12. With Continuous Operation:

- a. Perform a zero check (item 8.d) hourly
- b. Test a standard or optionally a control/sample (item 9) in the 600K to 800K range hourly in triplicate and determine the average, must be within 5% of the original established instrument average value (optionally, within 10% of original DMSCC average)
- c. Maintain records

### 13. Routine maintenance

- a. Maintain records

## REPORTS

### 14. Computing and Reporting Counts

- a. Count obtained x 1000 is the cell count/mL milk \_\_\_\_\_
- b. In reporting electronic somatic cell counts (ESCC/mL); record only first two left hand digits, raising second digit to next higher number when third digit is 6 or more \_\_\_\_\_
- c. Report the two left hand digits (rounded) \_\_\_\_\_
  - 1. If the third digit is 5 the second digit is rounded by the following rule \_\_\_\_\_
    - a. When the second digit is odd round up, raise the second digit by 1 (odd up, 235 to 240) \_\_\_\_\_
    - b. When the second digit is even round down, delete the 5 and report the second digit as is (even down, 225 to 220) \_\_\_\_\_
- d. If count on instrument is < 100 report as < 100,000 ESCC/mL \_\_\_\_\_
- e. If goat milk is over the regulatory limit, follow confirmation procedure in PMO \_\_\_\_\_