

Facility/Laboratory Name: \_\_\_\_\_ **INCUBATOR TEMPERATURE LOG**

**OFFICIAL RECORDS FOR THE PERIOD** \_\_\_\_\_ **TO** \_\_\_\_\_

Make/ Model: \_\_\_\_\_ Serial Number: \_\_\_\_\_ Unit Id#/ Identification#: \_\_\_\_\_ Temperature Range Used : \_\_\_\_\_

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
AM Top																																
AM Bottom																																
PM Top																																
PM Bottom																																

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
AM Top																																
AM Bottom																																
PM Top																																
PM Bottom																																

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
AM Top																																
AM Bottom																																
PM Top																																
PM Bottom																																

1. Size adequate for work load
2. Maintains appropriate temperature range for plate incubation.
3. Not to be used to store food or drink for consumption.
4. Record temperature (corrected) daily top and bottom shelves if applicable, in AM and PM, thermometer with bulb immersed in antifreeze liquid ( in sealed containers)
5. Record temperature daily to the nearest 0.1°C, maintain records. Note any corrective actions taken.

Example:

32.5°C MH
--------------