

Facility/Laboratory	Nama:
raciiity/Laboratory	Traffic.

STOCK PHOSPHATE BUFFER

Prepared on-site

Type: KH ₂ PO ₄ or MgCl ₂	Date Prepared	Date Autoclaved	Size of Storage Containers(mL)	Number of Storage Containers	Final pH	Analyst ID# or Initials

- 1. Final storage containers size less than or equal to 100ml.
- 2. Final pH 7.2 ± 0.2
- 3. Prepare dilution water with 1.25ml stock buffer/liter of MS water.
- 4. Label storage containers with date prepared and store after autoclaving in refrigerator.
- 5. Discard any buffer storage containers if discoloration occurred during storage period.
- 6. Do not use any buffer solution that may be questionable.