



Facility/Laboratory Name: \_\_\_\_\_

**MONTHLY MILK/MEDIA ELECTRONIC BALANCE CHECK RECORDS**

Year: \_\_\_\_\_

Make/ Mode/ Type: \_\_\_\_\_ Serial # or ID#: \_\_\_\_\_

Date(s) Serviced: \_\_\_\_\_

Date	Analysts ID# or Initials	Actual Scale Readings							Comments
		1 gm	5 gm	10 gm	25 gm	50 gm	100 gm	150 gm	

1. Electronic only, sensitive to 0.1g for general laboratory purposes and to 0.0001g for pipettor checks and antibiotics.
2. Checked monthly with Class S or S1, or equivalent ASTM 1, 2, or 3, weights. Certificate of authenticity or other verification, required.
3. Check weights that correspond to normal use of balance and maintain records. Record actual display reading of balance.
4. Checked at least annually by a qualified representative for good working order with proof of check in laboratory.