

Facility/Laboratory Name:			
	TEST KIT SUITABILITY	CHECK FOR DRUG RESIDUE TE	STING
Test Method Used:			

Date	Lot	Lot Expiration Da	Date	ate Date Start	Control Results/Interpretation			Analyst ID	
Received	Number	Date	Tested <sup>1</sup>	Using <sup>2</sup>	POSITIVE CONTROL		NEGATIVE CONTROL		or Initials
					Result	Suitability Check Date <sup>3</sup>	Result	Suitability Check Date <sup>3</sup>	
						Check Date <sup>3</sup>		Check Date <sup>3</sup>	
					Interpretation		Interpretation		

- 1. 'Date tested' is the date the first positive and negative control is run to check the suitability of the new lot.
- 2. 'Date Start Using' is the date this lot of kits is put into use, testing controls & trucks.
- 3. Suitability check date the POS/NEG control that was used to check the new kit was first tested and found suitable.