



Facility/Laboratory Name: _____

Microbic Density Records- Period _____ to _____ Year: _____

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
AM																																
PM																																

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

1. Microbic density of air less than (\leq) 15 colonies/plate or less than (\leq) 10 colonies/film in 15 min exposure.
2. The air control plate must be the first plate poured immediately **before** samples are shaken and must be located in the area of the plating activity (not off to the side).
3. For Petrifilm™ Aerobic Count Method: The film must be the first one prepared **before** samples are shaken and must be located in the area of activity (not off to the side).